

# FaithScriptor Reverse Transcriptase

Catalog No.: 2202-050 Size: 50 rxns

## Storage Condition

-20°C

## Components

Component	Amount
FaithScriptor Reverse Transcriptase(200 U/μl)	50 μl
5X RT Reaction Buffer	500 μl
0.1 M DTT	100 μl

## Description

FaithScriptor Reverse Transcriptase is a third-generation M-MLV reverse transcriptase obtained through genetic modification and recombinant technology. This enzyme has removed RNase H activity relative to wild-type M-MLV reverse transcriptase and significantly improved thermal stability (optimal reaction temperature is 50°C), thereby greatly enhancing the specificity and length of first-strand cDNA synthesis (up to 12 kb), and improving tolerance to complex RNA secondary structures.

## Definition of Activity Unit

The enzyme quantity needed to incorporate 1 nmol of [<sup>3</sup>H] dTTP in 10 minutes at 37°C using Poly(A)-Oligo(dT) as template/primer is defined as 1 unit of activity.

## Quality Control Assays

### Residual Host DNA

The product was tested by TaqMan qPCR with primers specific for the *E.coli* 16S rDNA, and the results show that the *E.coli* genome residues less than 1 copy/10 copies.

### Endonuclease Activity

The product was tested in a reaction containing a supercoiled plasmid DNA substrate. After incubation for 4 hours at 37°C, there was no significant change of the DNA substrate by agarose gel electrophoresis.

### Exonuclease Activity

The product was tested in a reaction containing DNA substrate. After incubation for 16 hours at 37°C, there was no significant change of the DNA substrate by agarose gel electrophoresis.

## Protocol

### First-strand cDNA synthesis:

1. Prepare the following reaction system on ice:

Reagent	Amount
Primer	X μl
Oligo(dT) <sub>20</sub>	The final concentration is 2.5 μM
Or Random Primer	The final concentration is 2.5 ng/μl
Or Gene-specific Primers	The final concentration is 0.25 μM
Template RNA <sup>a</sup>	50 ng~1 μg/20 μl
5X RT Reaction Buffer	4 μl
0.1 M DTT	1 μl
FaithScriptor Reverse Transcriptase (200 U/μl)	1 μl
dNTP Mix (10 mM Each)	1 μl
(Optional) RNase Inhibitor (40 U/μl)	1 μl
Nuclease-Free Water	To 20 μl

a. It is recommend to use high-quality RNA extracted using a kit that removes genomic DNA contamination as a template.

2. Mix gently and spin down.

3. If using Oligo(dT)<sub>20</sub> or gene-specific primers, incubate at 50°C for 30 min; If using random primers, first incubate at 25°C for 5 min, followed by incubation at 50°C for 30 minutes.

✦ Note: If the desired cDNA is less than 3 kb, the incubation time can be shortened to 15 minutes.

4. Terminate the reaction by incubating at 85°C for 5 minutes.

5. Place the cDNA solution on ice for use in subsequent experiments.

✦ Note: The cDNA solution can be stored at -20°C for up to six months, Long-term storage is recommended at -80°C.