

Serological Pipettes

Catalog No.:	1002-500	1003-500	1004-250
Pipet Volume:	5 mL	10 mL	25 mL
Size:	500 pcs/case	500 pcs/case	250 pcs/case

Purpose: Use for liquid transfer

Material:

- **Pipette:** GPPS , Compliance with USP Class VI
- **Pipette:** PO(polyolefin)
- **Color:** Clear




Features:

- Color-coded pipettes and packaging for convenient volume selection
- Individually packed pipettes in either paper/plastic or plastic peel for high sterility assurance
- High contrast graduations for easy-to-read volumes
- Ensures dispensing accuracy +/-2%
- Reverse graduations aid in determining dispensing volume
- Prevent aerosol or liquid contamination of the pipetting device while further reducing sample-to-sample contamination risks
- **Material:** Ultra clear, 100% virgin polystyrene for maximum clarity
- **Contaminant Free:** RNase/DNase free, Human DNA free, Non-Cytotoxic, Non-Pyrogenic
- **Sterility:** E-Beam irradiation to SAL of 10⁻⁶
- **Shelf Life:** 3 years from Manufacturing Date

DISPOSABLE SEROLOGICAL PIPETTE

Pipette products series:

BioCluster pipettes are made of medical grade exclusively high-grade polystyrene(GPPS) which are excellent for clear observation and reducing liquid attachment on the pipette surface to assure accurate delivery. The universal design makes the pipettes suitable for most brands of pipette-aids. The quality is guaranteed by the control and management of the production processing according to the ISO9001 and ISO13485. Six different sizes sterilized or non sterilized are provided: 1, 2, 5, 10, 25,50ml. Graduations are calibrated for accurate dispensing to within $\pm 2\%$.



50mL

25mL

10mL

5mL

2mL

1mL